

Technical Data Sheet Hypodermic Syringes



Sanah Luer Slip and Luer Lock hypodermic syringes deliver exceptionally smooth sliding and injection performance thanks to their precise plunger action, ensuring accurate measurements for a wide range of procedures.

Product name	Hypodermic syringes		
Characteristics	Sterile Latex-free Non-pyrogenic & non-toxic Individual soft blister packs Clear barrel with precise scale markings Compatible with standard hypodermic needles Low-dose sharp-tip stopper included in 1 mL syringes		
Material	Polypropylene plastic (PP)		
Quality compliance standards	ISO 7886-1	ISO 7864	ISO 80369-7
	Manufactured in an ISO 13485 and MDSAP certified facility		
Shelf life	5 years		
Storage conditions	Store in a dry place at room temperature		
Use	For veterinary use only		
Manufacturing origin	Made in China		

Supplier code	Description	Tip	Format
SA-SYRLS-1	Hypodermic syringes LS 1 mL	Luer Slip	100 units
SA-SYRLS-3	Hypodermic syringes LS 3 mL	Luer Slip	100 units
SA-SYRLS-5	Hypodermic syringes LS 5 mL	Luer Slip	100 units
SA-SYRLS-10	Hypodermic syringes LS 10 mL	Luer Slip	100 units
SA-SYRLL-1	Hypodermic syringes LL 1 mL	Luer Lock	100 units
SA-SYRLL-3	Hypodermic syringes LL 3 mL	Luer Lock	100 units
SA-SYRLL-5	Hypodermic syringes LL 5 mL	Luer Lock	100 units
SA-SYRLL-10	Hypodermic syringes LL 10 mL	Luer Lock	100 units
SA-SYRLL-20	Hypodermic syringes LL 20 mL	Luer Lock	40 units
SA-SYRLL-30	Hypodermic syringes LL 30 mL	Luer Lock	40 units
SA-SYRLL-60	Hypodermic syringes LL 60 mL	Luer Lock	30 units