

Technical Data Sheet

Woven Gauze Sponges



Sanah woven gauze sponges are premium high-density sponges made from 100% cotton, providing high absorbency. They are ideal for daily veterinary activities requiring a more heavy-duty and durable sponge resistant to tears. Perfect for general cleaning, prepping intact skin, disinfection, wound debriding, and applying pressure.

Product name	Woven non-sterile gauze sponges
Characteristics	Non-sterile Latex-free Single use
Number of layers	4 ply
Material	100% cotton
Thread Count	19 x 15 mesh
Absorbency	Retention: ≥ 9 g/g
	Submersion: ≤ 10 sec
Quality compliance standards	EN 14079-2003 Manufactured in an ISO 13485 certified facility
Shelf life	5 years
Storage conditions	Store in a dry place at room temperature
Use	For veterinary use only
Manufacturing origin	Made in China

Supplier code	Description	Format
SA-GWN4P3X3	Gauze sponges woven non-sterile 4 ply 3 x 3 in.	Pack of 100