

Technical Data Sheet

Vinyl Veterinary Examination Gloves



Sanah vinyl examination gloves are high-quality products that are cost-effective when executing non-medical veterinary tasks requiring a low level of protection.

Product name	Vinyl veterinary examination gloves
Characteristics	Non-sterile Latex-free Powder-free Ambidextrous
Material	Vinyl 3.5 mil (clear)
Thickness	Finger ≥ 0.08 mm
	Palm ≥ 0.08 mm

Dimensions			
Size	Length (mm)	Width (mm)	Average Weight (g)
S	≥ 230	85 ± 5	4.6
M	≥ 230	95 ± 5	5.0
L	≥ 230	105 ± 5	5.5
XL	≥ 230	115 ± 5	6.0

Physical property characteristics		
Tensile strength	Before aging	After aging
	≥ 11 MPa	≥ 11 MPa
Ultimate elongation	Before aging	After aging
	$\geq 300\%$	$\geq 300\%$

Powder residue	≤ 2.0 mg per glove
-----------------------	-------------------------

Technical Data Sheet

Vinyl Veterinary Examination Gloves



Quality compliance standards	ASTM F1671 ASTM D5250 EN 455 EN 374 Manufactured in an ISO 13485 and MDSAP certified facility
Shelf life	5 years
Storage conditions	Store in a dry place at room temperature
Use	For veterinary use only
Manufacturing origin	Made in China

Supplier code	Description	Format
SA-GLVNV-S	Vinyl veterinary examination gloves small (S)	100 gloves
SA-GLVNV-M	Vinyl veterinary examination gloves medium (M)	100 gloves
SA-GLVNV-L	Vinyl veterinary examination gloves large (L)	100 gloves
SA-GLVNV-XL	Vinyl veterinary examination gloves extra large (XL)	100 gloves