

Technical Data Sheet

Non-Woven Gauze Sponges



Sanah non-woven gauze sponges are premium products made from 70% viscose and 30% polyester. The high thread density of these 30 gsm sponges offers high absorbency, making them ideal for daily veterinary activities. The high density also means that less layers and less sponges are needed to accomplish your tasks.

Product name	Non-woven non-sterile gauze sponges
Characteristics	Non-sterile Latex-free Single use
Number of layers	4 ply
Material	70% viscose + 30% polyester (30 gsm)
Absorbency	Retention: ≥ 10 g/g
	Submersion: ≤ 10 sec
pH	$6 \leq \text{pH} \leq 8$ (neutral)
Tensile strength	MD $\geq 33\text{N}$ CD $\geq 6\text{N}$
Quality compliance standards	EN 14079-2003 Manufactured in an ISO 13485 certified facility
Shelf life	5 years
Storage conditions	Store in a dry place at room temperature
Use	For veterinary use only
Manufacturing origin	Made in China

Supplier code	Description	Format
SA-GNN4P2X2	Gauze sponges non-woven non-sterile 4 ply 2 x 2 in.	Pack of 200
SA-GNN4P3X3	Gauze sponges non-woven non-sterile 4 ply 3 x 3 in.	Pack of 200
SA-GNN4P4X4	Gauze sponges non-woven non-sterile 4 ply 4 x 4 in.	Pack of 200