

Technical Data Sheet

Hypodermic Syringes with Needles



Sanah hypodermic syringes with needles provide a simple yet high-quality solution for everyday veterinary procedures. The needles are manufactured with ultra-sharp tri-beveled tips for smooth and precise penetration, minimizing animal discomfort. They are securely pre-fixed to their syringe, streamlining workflow and improving operational efficiency.

Product name	Hypodermic syringes with needles
Overall characteristics	<ul style="list-style-type: none"> Sterile Latex-free Non-pyrogenic Non-toxic Individual soft blister packs
Syringe characteristics	<ul style="list-style-type: none"> Clear barrel with precise scale markings Low-dose sharp-tip stopper included in 0.5 mL and 1 mL syringes Made of polypropylene plastic (PP)
Needle characteristics	<ul style="list-style-type: none"> Tri-beveled edge and ultra-sharp tip minimizing patient discomfort Made of stainless steel Colour-coded for easy size identification
Quality compliance standards	<ul style="list-style-type: none"> ISO 7864 ISO 80369-7 ISO 6009 ISO 9626 Manufactured in an ISO 13485 and MDSAP certified facility
Shelf life	5 years
Storage conditions	Store in a dry place at room temperature
Security and disposal	<ul style="list-style-type: none"> Do not resheild used needles Dispose of sharps after single use according to facility protocol to avoid the risk of infection
Use	For veterinary use only
Manufacturing origin	Made in China

Technical Data Sheet

Hypodermic Syringes with Needles



Supplier code	Description	Tip	Colour	Format
SA-SYRNLS-05M2712	Hypodermic syringes LS 0.5 mL with 27G x ½" (0.4 x 13 mm) needles	Luer Slip		100 units
SA-SYRNLS-1M2712	Hypodermic syringes LS 1 mL with 27G x ½" (0.4 x 13 mm) needles	Luer Slip		100 units
SA-SYRNLS-1M2558	Hypodermic syringes LS 1 mL with 25G x ⅝" (0.5 x 16 mm) needles	Luer Slip		100 units
SA-SYRNLS-3M221	Hypodermic syringes LS 3 mL with 22G x 1" (0.7 x 25 mm) needles	Luer Slip		100 units
SA-SYRNLS-3M2234	Hypodermic syringes LS 3 mL with 22G x ¾" (0.7 x 19 mm) needles	Luer Slip		100 units
SA-SYRNLL-1M2558	Hypodermic syringes LL 1 mL with 25G x ⅝" (0.5 x 16 mm) needles	Luer Lock		100 units
SA-SYRNLL-3M2558	Hypodermic syringes LL 3 mL with 25G x ⅝" (0.5 x 16 mm) needles	Luer Lock		100 units
SA-SYRNLL-3M231	Hypodermic syringes LL 3 mL with 23G x 1" (0.6 x 25 mm) needles	Luer Lock		100 units
SA-SYRNLL-3M221	Hypodermic syringes LL 3 mL with 22G x 1" (0.7 x 25 mm) needles	Luer Lock		100 units
SA-SYRNLL-3M2234	Hypodermic syringes LL 3 mL with 22G x ¾" (0.7 x 19 mm) needles	Luer Lock		100 units
SA-SYRNLL-3M211	Hypodermic syringes LL 3 mL with 21G x 1" (0.8 x 25 mm) needles	Luer Lock		100 units
SA-SYRNLL-3M20112	Hypodermic syringes LL 3 mL with 20G x 1 ½" (0.9 x 38 mm) needles	Luer Lock		100 units
SA-SYRNLL-3M201	Hypodermic syringes LL 3 mL with 20G x 1" (0.9 x 25 mm) needles	Luer Lock		100 units