

Technical Data Sheet Insulin Syringes



Sanah insulin syringes are built to meet the demands of busy veterinary practices. They are available in two calibrations, ensuring accurate dosing for diabetic animals. They are offered in individual blister packs suitable for in-clinic use as well as client home use.

Product name	Insulin syringes
Overall characteristics	<ul style="list-style-type: none"> Sterile Latex-free Non-pyrogenic Non-toxic Individual soft blister packs
Syringe characteristics	<ul style="list-style-type: none"> Bold and clear graduation for accurate dosage Half-unit markings on 0.3 mL and 0.5 mL syringes Smooth plunger action for greater filling accuracy Made of polypropylene plastic (PP)
Needle characteristics	<ul style="list-style-type: none"> Tri-beveled edge and ultra-sharp tip minimizing patient discomfort Made of stainless steel
Quality compliance standards	<ul style="list-style-type: none"> ISO 9626 ISO 10993-7 Manufactured in an ISO 13485 and MDSAP certified facility
Shelf life	5 years
Storage conditions	Store in a dry place at room temperature
Security and disposal	<ul style="list-style-type: none"> Do not reshield used needles Dispose of sharps after single use according to facility protocol to avoid the risk of infection
Use	For veterinary use only
Manufacturing origin	Made in China

Technical Data Sheet

Insulin Syringes



Supplier code	Description	Colour	Format
SA-SYRIN40-03M2912	Insulin syringes U-40 0.3 mL with 29G x 1/2" (0.33 x 13 mm) needles	Red	100 units
SA-SYRIN40-05M2912	Insulin syringes U-40 0.5 mL with 29G x 1/2" (0.33 x 13 mm) needles	Red	100 units
SA-SYRIN40-1M2912	Insulin syringes U-40 1 mL with 29G x 1/2" (0.33 x 13 mm) needles	Red	100 units
SA-SYRIN40-1M2812	Insulin syringe U-40 1 mL with 28G x 1/2" (0.36 x 13 mm) needles	Teal	100 units
SA-SYRIN100-03M31516	Insulin syringes U-100 0.3 mL with 31G x 5/16" (0.25 x 8 mm) needles	White	100 units
SA-SYRIN100-03M2912	Insulin syringes U-100 0.3 mL with 29G x 1/2" (0.33 x 13 mm) needles	Red	100 units
SA-SYRIN100-05M2912	Insulin syringes U-100 0.5 mL with 29G x 1/2" (0.33 x 13 mm) needles	Red	100 units
SA-SYRIN100-05M2812	Insulin syringes U-100 0.5 mL with 28G x 1/2" (0.36 x 13 mm) needles	Teal	100 units