

Technical Data Sheet Hypodermic Needles



Sanah needles are designed for use with animals in routine veterinary procedures. They are precision-manufactured for consistent sharpness and smooth penetration, promoting controlled injections and reduced tissue trauma.

Product name	Hypodermic needles
Characteristics	<ul style="list-style-type: none"> Sterile Latex-free Non-pyrogenic Non-toxic Individual soft blister packs Compatible with standard hypodermic syringes Tri-beveled edge and ultra-sharp tip minimizing patient discomfort Colour-coded for easy size identification
Material	Stainless steel
Quality compliance standards	<ul style="list-style-type: none"> ISO 7884 ISO 6009 ISO 9626 ISO 80369-7 Manufactured in an ISO 13485 and MDSAP certified facility
Shelf life	5 years
Storage conditions	Store in a dry place at room temperature
Security and disposal	<ul style="list-style-type: none"> Do not reshield used needles Dispose of sharps after single use according to facility protocol to avoid the risk of infection
Use	For veterinary use only
Manufacturing origin	Made in China

Technical Data Sheet

Hypodermic Needles



Supplier code	Description	Colour	Format
SA-NEED-161	Hypodermic needles 16G x 1" (1.6 mm x 25 mm)		100 units
SA-NEED-16112	Hypodermic needles 16G x 1 ½" (1.6 mm x 38 mm)		100 units
SA-NEED-181	Hypodermic needles 18G x 1" (1.2 mm x 25 mm)		100 units
SA-NEED-18112	Hypodermic needles 18G x 1 ½" (1.2 mm x 38 mm)		100 units
SA-NEED-201	Hypodermic needles 20G x 1" (0.9 mm x 25 mm)		100 units
SA-NEED-20112	Hypodermic needles 20G x 1 ½" (0.9 mm x 38 mm)		100 units
SA-NEED-211	Hypodermic needles 21G x 1" (0.8 mm x 25 mm)		100 units
SA-NEED-2234	Hypodermic needles 22G x ¾" (0.7 mm x 19 mm)		100 units
SA-NEED-221	Hypodermic needles 22G x 1" (0.7 mm x 25 mm)		100 units
SA-NEED-22112	Hypodermic needles 22G x 1 ½" (0.7 mm x 38 mm)		100 units
SA-NEED-2334	Hypodermic needles 23G x ¾" (0.6 mm x 19 mm)		100 units
SA-NEED-231	Hypodermic needles 23G x 1" (0.6 mm x 25 mm)		100 units
SA-NEED-2558	Hypodermic needles 25G x 5/8" (0.5 mm x 16 mm)		100 units
SA-NEED-2612	Hypodermic needles 26G x ½" (0.45 mm x 13 mm)		100 units
SA-NEED-2712	Hypodermic needles 27G x ½" (0.4 mm x 13 mm)		100 units