

Technical Data Sheet

Oral Syringes with Caps



Sanah oral syringes are designed for the accurate administration of liquid medications, supplements, and nutritional support in animal patients. They are offered in sterile individual packages suitable for in-clinic use or client home use.

Product name	Oral syringes with caps
Characteristics	<ul style="list-style-type: none"> Sterile Latex-free Non-pyrogenic Non-toxic Individual soft blister packs Smooth clear barrel with precise markings and purple plunger Centric or catheter tip
Material	Polypropylene plastic (PP)
Quality compliance standards	<ul style="list-style-type: none"> ISO 7886-1 ISO 10993 Manufactured in an ISO 13485 and MDSAP certified facility
Shelf life	5 years
Storage conditions	Store in a dry place at room temperature
Use	For veterinary use only
Manufacturing origin	Made in China

Supplier code	Description	Tip	Format
SA-SYROR-1	Oral syringes with caps 1 mL	Centric	100 units
SA-SYROR-3	Oral syringes with caps 3 mL	Centric	100 units
SA-SYROR-5	Oral syringes with caps 5 mL	Catheter	100 units
SA-SYROR-10	Oral syringes with caps 10 mL	Catheter	48 units