

Technical Data Sheet

Latex Veterinary Examination Gloves



Sanah latex examination gloves are high-quality products offering maximum tactile sensitivity as well as stretch and strength. They are ideal for daily veterinary activities requiring precise or delicate handling.

Product name	Latex veterinary examination gloves		
Characteristics	Non-sterile Powder-free Ambidextrous Textured fingers for secure grip Contain latex which may cause allergic reactions		
Material	Natural latex rubber (ivory)		
Thickness	Finger	0.12 ± 0.02 mm	
	Palm	0.10 ± 0.02 mm	
Inner surface	Powder-free, chlorinated for easy donning		

Dimensions			
Size	Length (mm)	Width (mm)	Weight (g)
XS	≥ 230	≤ 80 ± 5	3.5 ± 0.3
S	≥ 230	80 ± 5	4.4 ± 0.3
M	≥ 230	95 ± 5	5.0 ± 0.3
L	≥ 230	110 ± 5	5.6 ± 0.3
XL	≥ 230	≥ 110	6.3 ± 0.3

Force at break (N)	Before aging	After aging
	≥ 6N (median)	≥ 6N (median)

Protein content max.	≤ 50 µg/dm ²
Powder residue	≤ 2.0 mg per glove

Technical Data Sheet

Latex Veterinary Examination Gloves



Quality compliance standards	ASTM F1671 ISO 11193-1 ISO 21420:2020 EN 455 EN 374 Manufactured in an ISO 13485 and MDSAP certified facility
Shelf life	5 years
Storage conditions	Store in a dry place at room temperature Keep out of direct sunlight and away from sources of ozone or ignition
Use	For veterinary use only
Manufacturing origin	Made in China

Supplier code	Description	Format
SA-GLVNL-XS	Latex veterinary examination gloves extra small (XS)	100 gloves
SA-GLVNL-S	Latex veterinary examination gloves small (S)	100 gloves
SA-GLVNL-M	Latex veterinary examination gloves medium (M)	100 gloves
SA-GLVNL-L	Latex veterinary examination gloves large (L)	100 gloves
SA-GLVNL-XL	Latex veterinary examination gloves extra large (XL)	100 gloves